



## Bloodless Castration Tool



Weiyou® Bloodless castration tools is made of high quality stainless steel. Weiyou® Bloodless castration tools is used for castration of bulls, calves, rams, pigs and colts. As the manufacturer in China, Weiyou® has mature production lines and supply chains, and always supply best quality product, best favourable price and best service.

## Product Description

Weiyou® Bloodless castration tools is made of high quality stainless steel. Weiyou® Bloodless castration tools is used for castration of bulls, calves, rams, pigs and colts. Weiyou® Bloodless castration tools is high degree of aesthetic and corrosion resistance, and be manufactured for optimal results & precision. OEM customization is acceptable. Free samples, and small quantity can be supplied.

## Specification

Item	Bloodless castration tools
Usage	For the humane and bloodless castration of livestock
Material	High quality stainless steel
MOQ	100pcs
Specification	37cm: For bulls 1-2 years old. 52cm: For large mature cattle.
Use For	Cattle, Horse Sheep, Goat, etc
Packing	Carton

## Features

- Non Slip Grip Premium Quality Handle
- Polish to High Standard Finish
- High Degree of Aesthetic and Corrosion Resistance

























