



Plastic-Steel Veterinary Syringes

Weiyou® Plastic-Steel Veterinary Syringes is made of high quality plastic (PC or TPX) with UV additive. As the manufacturer in China, Weiyou has mature production lines and supply chains, and always supply best quality product, best favourable price and best service.

Product Description

Weiyou® Plastic-Steel Veterinary Syringes is high quality plastic (PC or TPX) with UV additive. Plastic-Steel Veterinary Syringes is mainly used for Cattle, Horse Sheep, Goat and other animals. Various colors are available, and OEM customization is acceptable. Vaccination to prevent disease for all kinds of animals. Plastic-Steel Veterinary Syringes is continuous injection and accurate scale. Free samples, and small quantity can be supplied.

Specification

Item	Plastic-Steel Veterinary Syringes
Usage	Vaccination
Material	Plastic, Steel
MOQ	100pcs
Specification	5ml, 10ml, 20ml, 30ml, 50ml
Use For	Cattle, Horse Sheep, Goat
Advantage	Continuous injection, accurate scale
Color	Orange, Blue and other colors
Type	Vaccination to prevent disease
Packing	Carton

Features



Specification: 5ml, 10ml, 20ml, 30ml, 50ml
Durable syringe with luer lock adaptor;
Ergonomic design; Reusable;
Metal rod with engraved scale;
High quality plastic (PC or TPX) with UV additive;
Color of plastic piston, cover and handle are favorable;
Transparent barrel with engraved scale or transparent scale are favorable;
Sterilized by boiled water.








Ergonomically designed handle







Export type-50ml













Accurate scale and adjustable

